

Standard solution-

A solution of known concentration is called standard solution.

Titration;-

The process of adding the solution taken in the burette dropwise to a fixed volume of other solution taken in the conical flask till the endpoint ,is called titration.

End point-

The stage at which the reaction is just completed.

Indicators-

The endpoint is detected with the help of some substance is called indicators. They exhibit a marked colour change at the endpoint.

Acid -base titration

Two type;- i) **Acidimetry**-determination of amount of base in a solution by titrating with a standard acid solution.

ii) **Alkalimetry**;- determination of amount of acid in a acid solution by titrating against a standard base solution.

pH range of indicators-

orange red to golden yellow.

Indicator	pH Range	Colour	
		Acid Medium	Basic Medium
Phenolphthalene	8.3 - 10	Colourless	Pink
Methyl Orange	3.1 - 4.5	Pink	Yellow
Litmus	4.5 - 8.3	Red	Blue

Choice of indicators in acid -Base titration

- Titration of strong acid against strong base:-

Eg- $\text{HCl}/\text{HNO}_3/\text{H}_2\text{SO}_4$ VS NaOH/KOH

Methyl orange or Phenolphthalein

- Titration of weak Acid Against Strong Base

Eg- Oxalic acid/Acetic acid vs NaOH/KOH

Phenolphthalein

- **Titration of Strong acid against Weak Base**



Methyl orange

- **Titration of Weak acid against weak Base**



Titration has no sharp pH change ,none of the indicator cannot be used

Principle of Volumetric analysis;-

$$N_1V_1=N_2V_2$$

N_1 =normality of titrant

V_1 =volume of titrant

N_2 =normality of titrate

V_2 =volume of titrate

1. 500ml of a decinormal solution is diluted by adding 300ml of water. what is the normality of resulting solution?

$$V_1 = 24 \text{ ml}$$

$$N_1 = ?$$

$$V_2 = 20 \text{ ml}$$

$$N_2 = 0.1 \text{ N}$$

$$V_1 N_1 = V_2 N_2$$

$$\begin{aligned} \text{Then } N_1 &= V_2 N_2 / V_1 \\ &= 20 * 0.1 / 24 \\ &= 0.083 \text{ N} \end{aligned}$$

$$v = 400 \text{ ml}$$

$$E = 49 \text{ gm}$$

$$N = 0.083 \text{ N}$$

$$N = w * 1000 / E * V$$

$$\begin{aligned} \text{Then } w &= N * E * V / 1000 \\ &= 0.08 * 49 * 400 / 1000 \\ &= 1.63 \text{ gm} \end{aligned}$$

$$w=0.252\text{g}$$

$$v=200\text{ml}$$

$$E=45\text{g}$$

$$N=w*1000/E*v$$

$$=0.252*1000/45*200$$

$$=0.028\text{N}$$

Calculate the normality of KOH solution containing
1.4g KOH in 700ml?

$$w=1.4\text{g}$$

$$V=700\text{ml}$$

$$E=56\text{gm}$$

$$N=w1000/EV$$

$$=1,4*1000/56*700$$

$$=0.0357\text{N}$$

Calculate the pH of 0.02M NaOH?

$$[\text{OH}^-] = 0.02\text{M}$$

$$\begin{aligned}\text{pOH} &= -\log[\text{OH}^-] \\ &= -\log[0.02] \\ &= 1.69\end{aligned}$$

$$\text{pH} + \text{pOH} = \text{pK}_w = 14$$

$$\begin{aligned}\text{So pH} &= 14 - \text{pOH} \\ &= 14 - 1.69 \\ &= 12.30\end{aligned}$$